

Pen-type scintillation detector

Features:

- Built-in spectroscopy amplifier
- Built-in single channel analyzer
- Built-in high voltage generator
- Small in size

By using 10 mm diameter Pen-type photomultiplier tubes, small scintillation crystals (with or without an internal or external collimator) can provide information on activity, for example in node sentinel studies.

These instruments can be supplied complete with a sealed electronic unit which contains a high voltage generator (DC controlled), spectroscopic amplifier and a single channel analyzer.

Pen-type instruments provide as well an analog output for spectroscopic information (pulse height spectra), and a TTL (SCA) output signal selecting pulses in the region of interest.

Lower Level Discriminator (LLD) and Upper Level Discriminator (ULD) are provided with test points and the analog output can be monitored with any commercial ADC.

The detector assembly

The probe consists of a NaI(Tl) scintillation crystal 9.5 x 20 mm, optically coupled to a Pen-type photomultiplier tube (Hamamatsu R1635). The assembly is housed in an aluminum housing with a 1.5 m long output cable.



The signal and power supply junction box provides the high voltage to the probe and processes the analog signal from the PMT.

Electronic working principles

The photomultiplier tube is fed using a resistor chain with a total resistance value of 3.7 M Ω .

The high voltage is positive and the anode is AC coupled to the input of the amplifier via a DC decoupling capacitor. The high voltage is controlled by a DC input voltage which is taken from the positive 12 volt power supply.

The control voltage can be monitored at the test point on the side of the junction box. The reference voltage is internally regulated. The high voltage can be calculated using:

$$\text{HV[kV]} = \text{Control Voltage [V]} \times 2.50 \text{ kV}$$

The maximum high voltage is -1250 V and the unit is protected against short-circuit or overload and polarity reversal.

The amplifiers

The amplifiers require a stable +12 and -12 V power supply voltage which are generated internally. The signal from the anode of the PMT is amplified with a shaping amplifier.

Specifications:

Type number

9.5 B 20 / 0.4 – E1 – X + Electronics

The probe

NaI(Tl) crystal: 9.5 mm x 20 mm
(different shapes/ crystals on request)
FWHM @ 662 keV < 8%
FWHM @ 122 keV < 17%
PMT type Hamamatsu R 1635
Aluminum housing
Output cable: RG174 1.5 m long
Connector: Lemo FFB 0S 403 CLAC32

Electronics

Amplifier

Shaping time: 0.63 μ s
Output pulse: Bipolar
1 μ s rise time
1 μ s fall time
Max. pulse height: +8 V

Single Channel Analyzer

Output pulse: TTL (+5 V)
Pulse length: 1.5 μ s
ULL adjustment: 20 turn pot
LLD adjustment: 20 turn pot
Factory Calibration:
ULD: Maximum
LLD: 20 keV

Junction box connections

Signal in: Lemo ERA 0S 403 CLL
Signal out: :Lemo ERA 00 250 CTL
TTL out: Lemo ERA 00 250 CTL
Mating cabling connectors provided. No cable provided. Cables can be supplied on special request.

Power Requirements

Voltage: +12 V
Power consumption: < 300 mW
Connector: Lemo connector with 3 m cable with flying leads (Brown = +12 V, White is ground)

Overall Size: 112 mm x 31 mm x 62 mm

Scionix Holland B.V.
P.O. Box 143
3980 CC Bunnik
The Netherlands

Tel: + 31 30 657 03 12
Fax: +31 30 656 75 63
Website: www.scionix.nl
E-mail: sales@scionix.nl